

Characterization of *Anaeromyxobacter dehalogenans* and Related Metal-Reducing Bacteria by MALDI-TOF MS


Kerry E. Preston^{1,2}, John R. Barr¹, Adrian Woolfitt¹, Hercules Moura¹, Sara Henry², Qingzhong Wu², Youlboong Sung², Frank E. Loeffler²

¹Centers for Disease Control and Prevention, 4770 Buford Hwy, MS F47, Atlanta, GA 30341

²Georgia Institute of Technology, School of Civil and Environmental Engineering, 311 Ferst Drive, ES&T Bldg., Atlanta, GA 30332







Procedure



1. Sample preparation:
 - a. collect and centrifuge cells
 - b. wash cell pellet twice with sucrose buffer



2. Spot sample on MALDI plate



ANAEROMYXOBACTER




Classification

STANDARDS										Total # Spectra
		2CP-C (Fumarate)	2CP-C (Nitrate)	2CP-C (2-CP)	2CP-3	2CP-1	R	DS str BB1	Geo str SZ	
INPUT CLASSES	Class #	1	2	3	4	5	6	7	8	
2CP-C (Fumarate)	1									9
2CP-C (Nitrate)	2									9
2CP-C (2-CP)	3									9
2CP-3	4									9
2CP-1	5									9
R	6									9
DS str BB1	7									9
Geo str. SZ	8									9


Scale: 0% 100% (10% steps)

Classification with the Random Forest Algorithm: Nine processed spectra from each class were grouped successfully, manifested by the dark areas of this grid view. The darkness of areas in the grid is directly proportional to the strength of the correlation between each spectrum and the group into which it is classified by Random Forest. The grey areas represent similarity between different classes.

MALDI-TOF Mass Spectrometry



Spectra of Metal-Reducing Bacteria and *E. Coli*




Genus Specificity: Typical spectra of *Geobacter*, *Desulfuromonas*, and *Anaeromyxobacter*. These fingerprints are easily, visually distinguished. *E. Coli* is included as a Gram-negative control.

Organisms Analyzed


(with abbreviations)

Organism	Electron Donor	Electron Acceptor
<i>Anaeromyxobacter dehalogenans</i> str. 2CP-C 2CP-C (Fumarate)	Acetate	Fumarate
<i>Anaeromyxobacter dehalogenans</i> str. 2CP-C 2CP-C (Nitrate)	Acetate	Nitrate
<i>Anaeromyxobacter dehalogenans</i> str. 2CP-C 2CP-C (2-CP)	Acetate	2-Chlorophenol
<i>Anaeromyxobacter dehalogenans</i> str. 2CP-1 2CP-1	Acetate	Fumarate
<i>Anaeromyxobacter dehalogenans</i> str. 2CP-3 2CP-3	Acetate	Fumarate
<i>Anaeromyxobacter dehalogenans</i> str. R R	Acetate	Fumarate
<i>Desulfuromonas michiganans</i> str. BB1 DS str. BB1	Acetate	Fumarate
<i>Geobacter sp.</i> str. SZ Geo str. SZ	Acetate	Fumarate




Phylogeny: The metal reducers analyzed are closely genetically related.

Gel-view Representation: Spectra of 2CP-C Grown with Various Electron Acceptors



Growth Conditions: This gel-view translation facilitates visualization of the minor differences between mass spectra of 2CP-C grown with different substrates.



Class #	Class Name	Class #	Class Name
1	2CP-C (Fumarate)	5	2CP-1
2	2CP-C (Nitrate)	6	R
3	2CP-C (2-Chlorophenol)	7	DS str. BB1
4	2CP-3	8	Geo str. SZ

Alternative Representation of Random Forest Classification Results: 3-Dimensional view of the grouping of spectra, and association to the grid representation. The groups with greater similarity are found closer together in space in the figure. Each dot represents one spectra (9 spectra in each of 8 classes; one color for each class).

Conclusions

MALDI-TOF MS offers a rapid and simple alternative for fingerprint analysis of environmentally relevant microorganisms.

Mass spectral profiles can discriminate microorganisms by group and strain, providing the specificity required for environmental monitoring.

Spectra may offer information on the metabolic status of a particular strain at the time of sampling.

Future work will focus on:

- Identification of proteins; linking proteomics and genomics approaches
- Analysis of consortia and environmental samples

Acknowledgements

Department of Energy • Environmental Remediation Sciences Program of the Department of Energy
Natural and Accelerated Bioremediation Research, a program of the Department of Energy